



- High precision system
- Highly flexible: For HbA1c, CRP, Glucose and Hemoglobin
- NGSP/IFCC certified HbA1c test
- Up to four tests from one sample
- Superior results in laboratory quality
- Individual hematocrit correction
- Easy to use
- Colored touch screen
- Fully automated measurement

InnovaStar®

Efficient operation combined with convenient use

- Fully automated system – no manual steps required during measurement
- Single cartridge containing all reagents needed for testing
- Precalibrated tests – no time-consuming calibration – cost saving
- Measurement of up to four tests from one sample
- Software available in several languages (ENG, GER, FRA, PORT, RUSS, NOR)
- Additional languages available on request
- Intuitive software – no special knowledge needed
- Selectable units – no manual calculation
- Expandable parameter portfolio

Technical specifications	
Measurement principle	Photometric/Turbidimetric
Barcode identification	Automatic barcode scan for reagents; Barcode scanner for samples optionally
Thermostatic control	Measurement constantly at 37°C
Connectivity	RS232 connection; USB connection via adapter
Selectable units	e. g. for HbA1c mmol/mol, % DCCT
Dimensions	200 mm (W) x 150 mm (H) x 170 mm (D)
Weight	4 kg
User interface	2,8" colored touch screen

Handed over by:



Climate-neutral print
(Carbon neutral)
ID: DE 626-2XSFE6D
Printed on FSC-certified paper



DiaSys
Diagnostic Systems GmbH
Alte Strasse 9
65558 Holzheim
Germany

Phone: +49 6432 9146-0
Fax: +49 6432 9146-32
E-Mail: systems@diasys.de
www.diasys-diagnostics.com



8201161 March 2020

CHOOSING QUALITY.


InnovaStar®

POC Testing in Lab Quality

FASTER
ASSAY TIME!
HbA1c
< 5 min



Precise. Flexible. Easy.
The New Dimension for HbA1c, CRP,
Hemoglobin and Glucose POC Testing.



CHOOSING QUALITY.

InnovaStar® – POC Testing at a New Level


Easy to use, superior precision, rapid results: InnovaStar® brings POC testing to a new quality level. A fully automatic, compact analyzer for daily practice, providing the performance you would expect from a laboratory analyzer.



InnovaStar® – The Flexible Solution

Efficiency will be a decisive factor in the physicians’ ability to confront future challenges such as the rapid rise in the number of diabetics. This requires fast and user-friendly instruments that produce precise and reliable results. InnovaStar® offers a broad range of superior features, guaranteed to make your daily routine much easier.

InnovaStar® processes whole blood, serum and plasma, depending on the parameter used. One decisive advantage of InnovaStar® is its ability to carry out multiple tests from a single 10 µL whole blood sample.




Benefits of the 4 in 1 Principle:

- High flexibility – individual configuration of up to 4 tests from 1 sample
- Less harmful for the patient – blood withdrawal only once
- Time saving and cost-effective – immediate results, no multiple visits


Laboratory Quality in its Most Compact Form

InnovaStar® represents a system where reagent and device perfectly match. The underlying technology and reagents in established DiaSys quality are identical to the ones used on clinical laboratory analyzers. As for lab analyzers HbA1c determination on InnovaStar® is carried out according to DCCT/NGSP and IFCC guidelines. These preconditions ensure high quality results with reliable precision.

Method comparison
InnovaStar® vs HPLC



Method comparison
InnovaStar® vs lab analyzer



Precision in decision range				
	HbA1c IS	CRP IS	Glucose IS	Hemoglobin IS
Mean	6.6% DCCT	10.6 mg/L	125 mg/dL	9.12 g/dL
SD	0.096% DCCT	0.22 mg/L	1.59 mg/dL	0.06 g/dL
CV	1.44%	2.10%	1.27%	0.7%

Quick Guide

